



MASTERBOND

Proven and Reliable

MASTERBOND is a one-bottle 5th generation dental bonding system intended for durable bonding of composite materials to hard dental tissues and for surface preparation of prosthetic restorations prior to cementation.

Key Features

- **Durable Bond:**
MASTERBOND ensures strong and long-lasting bonding of composite fillings to hard dental tissues, translating into treatment stability and longevity.
- **Effective Reduction of Polymerization Shrinkage:**
MASTERBOND effectively reduces polymerization shrinkage, eliminating the risk of gaps and ensuring durable results.
- **Ideal Preparation of Prosthetic Surfaces:**
MASTERBOND prepares the surfaces of prosthetic restorations for reliable cementation, ensuring secure and aesthetic fixation of crowns, bridges, and other prosthetic works.
- **Single-Component and Universal:**
MASTERBOND is a single-component, universal bonding system that ensures ease of use and versatility in every dental procedure.
- **Convenient to Use:**
MASTERBOND features simple application, increasing working comfort and the final effect.

Indications for Use

MASTERBOND is used for bonding hard dental tissues with composite fillings, sealing the marginal gap, and reducing polymerization shrinkage. It is used for the preparation of prosthetic restoration surfaces prior to their cementation with dual-cure composite cements.



Volume

- Bottle of 4 ml

