

PRODUCT CATALOGUE



ARKONA

Created by Dentists for Dentists



FLOW-COLOR

A microhybrid flowable composite in 9 visible colours for special applications.

Available shades:
red, orange, yellow, green, blue, purple, pink, white, black



ETCHGEL

Dental etching system in gel consistency, easy to apply and remove. It contains 36% o-phosphoric acid.

Available in syringes:
2,6 g, 6,5 g, 13 g or 60 g





Available shades: A1, A2, A3, A3.5, OA2, OA3, B1, B2, C2, D2, D3, Transparent

CREATE

"Hard as nails"
micro-hybrid restorative
and light cure composite.
Designed for prosterior work.

6 g syringe

NEXT

"As soft as velvet"
modelling composite for dental
fillings,
micro-hybrid and light cure.
For anterior and prosterior work.

4 g syringe

blue  pink 



Available shades: A1, A2, A3, A3.5, OA2, OA3, B1, B2, C2, D2, D3, pink, blue, incisal, Multishade

FLOW-ART

User friendly, flowable,
micro-hybrid
restorative composite.

2 g syringes



Available in shades: A1, A2, A3, A3.5, OA2

SEALPRIM

Pit and fissure transparent
sealant with short (5-10 sec)
curing time. No adhesive is
required before.

Syringe 1 g








Available in 1 shade: Transparent

COLOURANT

A flowable composite for
characterisation and natural look
of restorative works.
Available in brownish shades,
gingival and pink.

Syringe 1 g



    
A B C GINGIVAL WHITE

COMCORD

The world's only dental fibre combined with composite in one. No flowable composite is required for application.

length - 60 mm + 2 x 5 mm
width $2,5 \pm 0,2$ mm
thickness $0,7 \pm 0,1$ mm



COMCORD SCISSORS

Designed for cutting the most durable dental fibres, with a special blade structure.

1 pcs



CHORDS

A versatile, light-curing dental fibre with multiple applications. Strong and flexible polyaramid bridge type fibres.

length: 10 cm
widths: 2 mm, 3 mm
or 4 mm



PERIO-FLUSH

A solution for rinsing periodontal pockets in a disposal syringe form. Domifene bromide-based formula.

Set of 10 x 3 ml syringe



GLAZE

A light-curing resin to increase the smoothness and gloss of composite and glass ionomer restorations. Transparent shade. The final step when working with fibres.

Bottle 3 ml



MASTERBOND

A 5th generation, 1-bottle, universal adhesive system.

Bottle 4 ml



Indications and related products:

CROWNS/ ENDOCROWNS

CREATE, COMCORD,
CHORDS, FLOW-ART,
COLOURANT, ETCHGEL,
MASTERBOND

SPLINTING, SPLINTING BRIDGES

COMCORD, CHORDS,
MASTERBOND, ETCHGEL,
CREATE, FLOW-ART,
COLOURANT, GLAZE

ROOT CANAL MARKING

FLOW-COLOR

BRIDGES, SPLINTING BRIDGES

CREATE, FLOW-ART,
COLOURANT,
COMCORD, CHORDS,
MASTERBOND,
ETCHGEL, GLAZE

ONLAYS/INLAYS

CREATE, FLOW-ART,
COLOURANT, COMCORD,
CHORDS, MASTERBOND,
ETCHGEL, GLAZE

TEMPORARY SPLINTING, RETAINERS

COMCORD, CHORDS,
FLOW-COLOR,
MASTERBOND,
ETCHGEL

BIG FILLING REINFORCEMENT

COMCORD,
CHORDS

CAVITY FILLINGS

NEXT, CREATE,
FLOW-ART, COLOURANT,
MASTERBOND,
ETCHGEL, GLAZE

BITE ADJUSTMENT

NEXT, CREATE,
FLOW-COLOR

PIT AND FISSURE SEALING

FLOW-COLOR,
SEALPRIM, ETCHGEL,
MASTERBOND

PERIODONTAL RINSING

PERIO-FLUSH,

Contact us:



ARKONA

Created by Dentists for Dentists

ARKONA Laboratorium Farmakologii Stomatologicznej

Nasutów 99C, 21-025 Niemce, Poland, EU

+48 664 407 296, +48 887 883 005

export@arkonadent.com



www.arkonadent.com/en

2026